

5'-AMP-activated protein kinase subunit beta-1

Catalog No: #AP80143



Package Size: #AP80143-1 50ug #AP80143-2 100ug #AP80143-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	5'-AMP-activated protein kinase subunit beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	AMPK subunit beta-1,5'-AMP-activated protein kinase 40 kDa subunit
Accession No.	Uniprot ID:P80387
Uniprot	P80387
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PRKAB1

Note: This product is for in vitro research use only