

# 5'-AMP-activated protein kinase catalytic subunit alpha-1

Catalog No: #AP80145

Package Size: #AP80145-1 50ug #AP80145-2 100ug #AP80145-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	5'-AMP-activated protein kinase catalytic subunit alpha-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AMPK subunit alpha-1,Acetyl-CoA carboxylase kinase,ACACA kinase,Hydroxymethylglutaryl-CoA reductase kinase,HMGCR kinase,Tau-protein kinase PRKAA1
Accession No.	Uniprot ID:Q5EG47
Uniprot	Q5EG47
GeneID	105787;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Prkaa1

Note: This product is for in vitro research use only