

Alpha-1A adrenergic receptor

Catalog No: #AP80151

Package Size: #AP80151-1 50ug #AP80151-2 100ug #AP80151-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-1A adrenergic receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Alpha-1A adrenoreceptor,Alpha-1A adrenoceptor,Alpha-1C adrenergic receptor,ADRA1C
Accession No.	Uniprot ID:P18130
Uniprot	P18130
GeneID	282134;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADRA1A

Note: This product is for in vitro research use only