

Beta-1 adrenergic receptor

Catalog No: #AP80153



Package Size: #AP80153-1 50ug #AP80153-2 100ug #AP80153-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Beta-1 adrenergic receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Beta-1 adrenoreceptor,Beta-1 adrenoceptor
Accession No.	Uniprot ID:Q28998
Uniprot	Q28998
GeneID	397355;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADRB1

Note: This product is for in vitro research use only