

Angiotensin-converting enzyme 2

Catalog No: #AP80176

Package Size: #AP80176-1 50ug #AP80176-2 100ug #AP80176-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin-converting enzyme 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ACE-related carboxypeptidase,Angiotensin-converting enzyme homolog,ACEH,Metalloprotease MPROT15,UNQ868/PRO1885
Accession No.	Uniprot ID:Q9BYF1
Uniprot	Q9BYF1
GeneID	59272;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ACE2

Note: This product is for in vitro research use only