

ATP-citrate synthase

Catalog No: #AP80223



Package Size: #AP80223-1 50ug #AP80223-2 100ug #AP80223-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ATP-citrate synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ATP-citrate (pro-S-)-lyase,Citrate cleavage enzyme
Accession No.	Uniprot ID:Q91V92
Uniprot	Q91V92
GeneID	104112;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Acly

Note: This product is for in vitro research use only