

RBP-Jkappa antibody

Catalog No: #22930



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RBP-Jkappa antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 258 and 498 of RBP-Jkappa
Target Name	RBP-Jkappa
Accession No.	Swiss-Prot:Q06330Gene ID:3516
Uniprot	Q06330
GeneID	3516;
Concentration	0.3mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 56kd
Western blotting: 1:500-1:3000
Immunofluorescence: 1:100-1:200


Images



Sample (30 ug of whole cell lysate)
A: Raji
7.5% SDS PAGE
Primary antibody diluted at 1: 1000



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using RBP-Jkappa antibody at 1: 100 dilution.

Note: This product is for in vitro research use only