

Aryl hydrocarbon receptor

Catalog No: #AP80367



Package Size: #AP80367-1 50ug #AP80367-2 100ug #AP80367-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aryl hydrocarbon receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ah receptor, Class E basic helix-loop-helix protein 76, bHLHe76, BHLHE76
Accession No.	Uniprot ID:P35869
Uniprot	P35869
GeneID	196;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: AHR

Note: This product is for in vitro research use only