

Alpha-actinin-3

Catalog No: #AP80448



Package Size: #AP80448-1 50ug #AP80448-2 100ug #AP80448-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Alpha-actinin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Alpha-actinin skeletal muscle isoform 3,F-actin cross-linking protein
Accession No.	Uniprot ID:Q08043
Uniprot	Q08043
GeneID	89;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ACTN3

Note: This product is for in vitro research use only