

Atrial natriuretic peptide receptor 2

Catalog No: #AP80492

Package Size: #AP80492-1 50ug #AP80492-2 100ug #AP80492-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Atrial natriuretic peptide receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Atrial natriuretic peptide receptor type B, ANP-B, Guanylate cyclase B, GC-B
Accession No.	Uniprot ID:Q6VWW5
Uniprot	Q6VWW5
GeneID	230103;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Npr2

Note: This product is for in vitro research use only