

Aldehyde dehydrogenase, dimeric NADP-preferring

Catalog No: #AP80499

Package Size: #AP80499-1 50ug #AP80499-2 100ug #AP80499-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aldehyde dehydrogenase, dimeric NADP-preferring
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Aldehyde dehydrogenase 4,Aldehyde dehydrogenase family 3 member A1,Dioxin-inducible aldehyde dehydrogenase 3,Ahd-4,Ahd4,Aldh3,Aldh4
Accession No.	Uniprot ID:P47739
Uniprot	P47739
GeneID	11670;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aldh3a1

Note: This product is for in vitro research use only