

Cytosol aminopeptidase

Catalog No: #AP80644

Package Size: #AP80644-1 50ug #AP80644-2 100ug #AP80644-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytosol aminopeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Leucine aminopeptidase 3,LAP-3,Leucyl aminopeptidase,Proline aminopeptidase,Prolyl aminopeptidase,Lapep
Accession No.	Uniprot ID:Q9CPY7
Uniprot	Q9CPY7
GeneID	66988;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lap3

Note: This product is for in vitro research use only