

Aquaporin-2

Catalog No: #AP80668



Package Size: #AP80668-1 50ug #AP80668-2 100ug #AP80668-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aquaporin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AQP-2,ADH water channel,Aquaporin-CD,AQP-CD,Collecting duct water channel protein,WCH-CD,Water channel protein for renal collecting duct
Accession No.	Uniprot ID:P56402
Uniprot	P56402
GeneID	11827;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aqp2

Note: This product is for in vitro research use only