

Apolipoprotein C-I

Catalog No: #AP80708



Package Size: #AP80708-1 50ug #AP80708-2 100ug #AP80708-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Apolipoprotein C-I
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Apo-CI,Apolipoprotein C1
Accession No.	Uniprot ID:P34928
Uniprot	P34928
GeneID	11812;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Apoc1

Note: This product is for in vitro research use only