

## AP-1 complex subunit mu-1

Catalog No: #AP80795

Package Size: #AP80795-1 50ug #AP80795-2 100ug #AP80795-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AP-1 complex subunit mu-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AP-mu chain family member mu1A,Adaptor protein complex AP-1 subunit mu-1,Adaptor-related protein complex 1 subunit mu-1,Clathrin assembly protein complex 1 mu-1 medium chain 1,Clathrin coat assembly protein AP47,Clathrin coat-associated protein AP47,Golgi adaptor HA1/AP1 adaptin mu-1 subunit,Mu-adaptin 1,Mu1A-adaptin,Cltnm
Accession No.	Uniprot ID:P35585
Uniprot	P35585
GeneID	11767;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ap1m1

Note: This product is for in vitro research use only