

Apolipoprotein C-II

Catalog No: #AP80809



Package Size: #AP80809-1 50ug #AP80809-2 100ug #AP80809-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Apolipoprotein C-II
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Apo-CII,Apolipoprotein C2
Accession No.	Uniprot ID:Q05020
Uniprot	Q05020
GeneID	105886299;11813;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Apoc2

Note: This product is for in vitro research use only