

Alpha-methylacyl-CoA racemase

Catalog No: #AP80843



Package Size: #AP80843-1 50ug #AP80843-2 100ug #AP80843-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-methylacyl-CoA racemase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	2-arylpropionyl-CoA epimerase,2-methylacyl-CoA racemase
Accession No.	Uniprot ID:P70473
Uniprot	P70473
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Amacr

Note: This product is for in vitro research use only