



A disintegrin and metalloproteinase with thrombospondin motifs 18

Catalog No: #AP80864

Package Size: #AP80864-1 50ug #AP80864-2 100ug #AP80864-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 18
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ADAM-TS 18,ADAMTS21
Accession No.	Uniprot ID:Q8TE60
Uniprot	Q8TE60
GeneID	170692;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADAMTS18

Note: This product is for in vitro research use only