

Axin-1

Catalog No: #AP80868



Package Size: #AP80868-1 50ug #AP80868-2 100ug #AP80868-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Axin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Axis inhibition protein 1,AXIN,AXN
Accession No.	Uniprot ID:O42400
Uniprot	O42400
GeneID	395786;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AXIN1

Note: This product is for in vitro research use only