

Calsequestrin-1

Catalog No: #AP81089

Package Size: #AP81089-1 50ug #AP81089-2 100ug #AP81089-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calsequestrin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Calsequestrin, skeletal muscle isoform
Accession No.	Uniprot ID:P31236
Uniprot	P31236
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CASQ1

Note: This product is for in vitro research use only