

Brain-specific angiogenesis inhibitor 3

Catalog No: #AP81156



Package Size: #AP81156-1 50ug #AP81156-2 100ug #AP81156-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Brain-specific angiogenesis inhibitor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	KIAA0550
Accession No.	Uniprot ID:O60242
Uniprot	O60242
GeneID	577;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BAI3

Note: This product is for in vitro research use only