

B-cell lymphoma 6 protein

Catalog No: #AP81169

Package Size: #AP81169-1 50ug #AP81169-2 100ug #AP81169-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B-cell lymphoma 6 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BCL-6,B-cell lymphoma 5 protein,BCL-5,Protein LAZ-3,Zinc finger and BTB domain-containing protein 27,Zinc finger protein 51,BCL5,LAZ3,ZBTB27,ZNF51
Accession No.	Uniprot ID:P41182
Uniprot	P41182
GeneID	604;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BCL6

Note: This product is for in vitro research use only