

# Calpain-3

Catalog No: #AP81233

Package Size: #AP81233-1 50ug #AP81233-2 100ug #AP81233-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                                                                                                    |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | Calpain-3                                                                                                                                          |
| Brief Description  | Recombinant Protein                                                                                                                                |
| Host Species       | E.coli                                                                                                                                             |
| Purification       | Greater than 90% by SDS-PAGE                                                                                                                       |
| Species Reactivity | Bovine                                                                                                                                             |
| Other Names        | Calcium-activated neutral proteinase 3,CANP 3,Calpain L3,Calpain p94,Muscle-specific calcium-activated neutral protease 3,New calpain 1,nCL-1,NCL1 |
| Accession No.      | Uniprot ID:P51186                                                                                                                                  |
| Uniprot            | P51186                                                                                                                                             |
| GeneID             | 281663;                                                                                                                                            |
| Tag Info           | His                                                                                                                                                |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                                                                                             |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CAPN3

Note: This product is for in vitro research use only