

Calcitonin gene-related peptide 1

Catalog No: #AP81270

Package Size: #AP81270-1 50ug #AP81270-2 100ug #AP81270-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcitonin gene-related peptide 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-type CGRP, Calcitonin gene-related peptide I, CGRP-I, Calc
Accession No.	Uniprot ID:Q99JA0
Uniprot	Q99JA0
GeneID	12310;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Calca

Note: This product is for in vitro research use only