

Brain-specific serine protease 4

Catalog No: #AP81276

Package Size: #AP81276-1 50ug #AP81276-2 100ug #AP81276-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Brain-specific serine protease 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BSSP-4, Serine protease 22, Serine protease 26, Tryptase epsilon, BSSP4, PRSS26, SP001LA, UNQ302/PRO343
Accession No.	Uniprot ID: Q9GZN4
Uniprot	Q9GZN4
GeneID	64063;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PRSS22

Note: This product is for in vitro research use only