

C-C motif chemokine 24

Catalog No: #AP81303



Package Size: #AP81303-1 50ug #AP81303-2 100ug #AP81303-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 24
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Eosinophil chemotactic protein 2,Eotaxin-2,Small-inducible cytokine A24
Accession No.	Uniprot ID:Q68Y86
Uniprot	Q68Y86
GeneID	445452;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCL24

Note: This product is for in vitro research use only