

Calreticulin

Catalog No: #AP81338



Package Size: #AP81338-1 50ug #AP81338-2 100ug #AP81338-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calreticulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CRP55,Calregulin,Endoplasmic reticulum resident protein 60,ERp60,HACBP,grp60,CRTC
Accession No.	Uniprot ID:P27797
Uniprot	P27797
GeneID	811;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CALR

Note: This product is for in vitro research use only