

Cyclin-dependent kinase inhibitor 1

Catalog No: #AP81438

Package Size: #AP81438-1 50ug #AP81438-2 100ug #AP81438-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase inhibitor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CDK-interacting protein 1,Melanoma differentiation-associated protein 6,MDA-6,p21,CAP20,CDKN1,CIP1,MDA6,PIC1,SDI1,WAF1
Accession No.	Uniprot ID:P38936
Uniprot	P38936
GeneID	1026;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDKN1A

Note: This product is for in vitro research use only