

# Cyclin-dependent kinase 4 inhibitor C

Catalog No: #AP81439

Package Size: #AP81439-1 50ug #AP81439-2 100ug #AP81439-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                        |
|--------------------|--------------------------------------------------------|
| Product Name       | Cyclin-dependent kinase 4 inhibitor C                  |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli                                                 |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Ms                                                     |
| Other Names        | Cyclin-dependent kinase 6 inhibitor,p18-INK4c,p18-INK6 |
| Accession No.      | Uniprot ID:Q60772                                      |
| Uniprot            | Q60772                                                 |
| GeneID             | 12580;                                                 |
| Tag Info           | His                                                    |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cdkn2c

Note: This product is for in vitro research use only