

T-cell-specific surface glycoprotein CD28

Catalog No: #AP81478

Package Size: #AP81478-1 50ug #AP81478-2 100ug #AP81478-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	T-cell-specific surface glycoprotein CD28
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CD28
Accession No.	Uniprot ID:P31042
Uniprot	P31042
GeneID	25660;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cd28

Note: This product is for in vitro research use only