

C-C motif chemokine 17

Catalog No: #AP81497

Package Size: #AP81497-1 50ug #AP81497-2 100ug #AP81497-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CC chemokine TARC, Small-inducible cytokine A17, Thymus and activation-regulated chemokine, SCYA17, TARC
Accession No.	Uniprot ID: Q92583
Uniprot	Q92583
GeneID	6361;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CCL17

Note: This product is for in vitro research use only