

eNOS Monoclonal Antibody(Mix)

Catalog No: #44212



Package Size: #44212-1 50ul #44212-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	eNOS Monoclonal Antibody(Mix)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Mix
Isotype	IgG1
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous eNOS protein.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein of eNOS
Target Name	eNOS
Other Names	NOS3; Nitric oxide synthase, endothelial; Constitutive NOS; cNOS; EC-NOS; Endothelial NOS; eNOS; NOS type III; NOSIII
Calculated MW	130kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of Rat Heart Tissue, diluted at 1:1000.



Western blot analysis of HUVEC Cell Lysate, diluted at 1:1000.

Wuhan Tongji Hospital

Note: This product is for in vitro research use only