

60 kDa heat shock protein, mitochondrial

Catalog No: #AP81512



Package Size: #AP81512-1 50ug #AP81512-2 100ug #AP81512-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	60 kDa heat shock protein, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	60 kDa chaperonin,Chaperonin 60,CPN60,Heat shock protein 60,HSP-60,Hsp60,HSP60
Accession No.	Uniprot ID:P49818
Uniprot	P49818
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSPD1

Note: This product is for in vitro research use only