

C-C chemokine receptor type 4

Catalog No: #AP81578

Package Size: #AP81578-1 50ug #AP81578-2 100ug #AP81578-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C chemokine receptor type 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C-C CKR-4,K5-5,CD194,CMKBR4
Accession No.	Uniprot ID:P51679
Uniprot	P51679
GeneID	1233;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCR4

Note: This product is for in vitro research use only