

Carboxypeptidase B2

Catalog No: #AP81595

Package Size: #AP81595-1 50ug #AP81595-2 100ug #AP81595-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carboxypeptidase B2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Carboxypeptidase R,CPR,Carboxypeptidase U,CPU,Thrombin-activable fibrinolysis inhibitor,TAFI,Tafi
Accession No.	Uniprot ID:Q9JHH6
Uniprot	Q9JHH6
GeneID	56373;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cpb2

Note: This product is for in vitro research use only