

Eotaxin

Catalog No: #AP81610



Package Size: #AP81610-1 50ug #AP81610-2 100ug #AP81610-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Eotaxin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C-C motif chemokine 11,Eosinophil chemotactic protein,Small-inducible cytokine A11,SCYA11
Accession No.	Uniprot ID:P51671
Uniprot	P51671
GeneID	6356;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCL11

Note: This product is for in vitro research use only