

T-cell surface glycoprotein CD8 alpha chain

Catalog No: #AP81621



Package Size: #AP81621-1 50ug #AP81621-2 100ug #AP81621-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	T-cell surface glycoprotein CD8 alpha chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CD8a
Accession No.	Uniprot ID:P31783
Uniprot	P31783
GeneID	281060;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD8A

Note: This product is for in vitro research use only