

C-C motif chemokine 4

Catalog No: #AP81639

Package Size: #AP81639-1 50ug #AP81639-2 100ug #AP81639-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Immune activation protein 2,ACT-2,Macrophage inflammatory protein 1-beta,MIP-1-beta,Protein H400,SIS-gamma,Small-inducible cytokine A4,Mip1b,Scya4
Accession No.	Uniprot ID:P14097
Uniprot	P14097
GeneID	20303;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ccl4

Note: This product is for in vitro research use only