

## Collagen alpha-1(XXVII) chain

Catalog No: #AP81663

Package Size: #AP81663-1 50ug #AP81663-2 100ug #AP81663-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                        |
|--------------------|--------------------------------------------------------|
| Product Name       | Collagen alpha-1(XXVII) chain                          |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli                                                 |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Rt                                                     |
| Accession No.      | Uniprot ID:Q80ZF0                                      |
| Uniprot            | Q80ZF0                                                 |
| GeneID             | 298101;                                                |
| Tag Info           | His                                                    |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Col27a1

Note: This product is for in vitro research use only