

Clusterin

Catalog No: #AP81746

Package Size: #AP81746-1 50ug #AP81746-2 100ug #AP81746-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Clusterin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Dimeric acid glycoprotein,DAG,Sulfated glycoprotein 2,SGP-2,Testosterone repressed prostate message 2,TRPM-2
Accession No.	Uniprot ID:P05371
Uniprot	P05371
GeneID	24854;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Clu

Note: This product is for in vitro research use only