

Calcium and integrin-binding protein 1

Catalog No: #AP81750

Package Size: #AP81750-1 50ug #AP81750-2 100ug #AP81750-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcium and integrin-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Calmyrin
Accession No.	Uniprot ID:Q17QE5
Uniprot	Q17QE5
GeneID	510141;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CIB1

Note: This product is for in vitro research use only