

Carbohydrate sulfotransferase 15

Catalog No: #AP81767

Package Size: #AP81767-1 50ug #AP81767-2 100ug #AP81767-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carbohydrate sulfotransferase 15
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	B-cell RAG-associated gene protein,N-acetylgalactosamine 4-sulfate 6-O-sulfotransferase,GalNAc4S-6ST,Brag,Galnac4s6st,Kiaa0598
Accession No.	Uniprot ID:Q91XQ5
Uniprot	Q91XQ5
GeneID	77590;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Chst15

Note: This product is for in vitro research use only