

Cell division cycle protein 23 homolog

Catalog No: #AP81769



Package Size: #AP81769-1 50ug #AP81769-2 100ug #AP81769-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Cell division cycle protein 23 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Anaphase-promoting complex subunit 8,APC8,Cyclosome subunit 8,ANAPC8
Accession No.	Uniprot ID:A1A4R8
Uniprot	A1A4R8
GeneID	512651;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CDC23

Note: This product is for in vitro research use only