

Cell growth regulator with EF hand domain protein 1

Catalog No: #AP81782

Package Size: #AP81782-1 50ug #AP81782-2 100ug #AP81782-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cell growth regulator with EF hand domain protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cell growth regulatory gene 11 protein,Hydrophobestin,Cgr11
Accession No.	Uniprot ID:Q8R1U2
Uniprot	Q8R1U2
GeneID	68567;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cgref1

Note: This product is for in vitro research use only