

CUB domain-containing protein 1

Catalog No: #AP81871

Package Size: #AP81871-1 50ug #AP81871-2 100ug #AP81871-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CUB domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Membrane glycoprotein gp140,Subtractive immunization M plus HEp3-associated 135 kDa protein,SIMA135,Transmembrane and associated with src kinases,CD318,TRASK,UNQ2486/PRO5773
Accession No.	Uniprot ID:Q9H5V8
Uniprot	Q9H5V8
GeneID	64866;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDCP1

Note: This product is for in vitro research use only