

Segment polarity protein dishevelled homolog DVL-1

Catalog No: #AP82089

Package Size: #AP82089-1 50ug #AP82089-2 100ug #AP82089-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Segment polarity protein dishevelled homolog DVL-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dishevelled-1,DSH homolog 1,Dvl
Accession No.	Uniprot ID:P51141
Uniprot	P51141
GeneID	13542;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dvl1

Note: This product is for in vitro research use only