

C-X-C motif chemokine 2

Catalog No: #AP82328

Package Size: #AP82328-1 50ug #AP82328-2 100ug #AP82328-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Macrophage inflammatory protein 2,MIP2,Mip-2,Mip2,Scyb2
Accession No.	Uniprot ID:P10889
Uniprot	P10889
GeneID	20310;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cxcl2

Note: This product is for in vitro research use only