

Coxsackievirus and adenovirus receptor homolog

Catalog No: #AP82349

Package Size: #AP82349-1 50ug #AP82349-2 100ug #AP82349-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Coxsackievirus and adenovirus receptor homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CAR,Car
Accession No.	Uniprot ID:P97792
Uniprot	P97792
GeneID	13052;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cxadr

Note: This product is for in vitro research use only