

CFTR (Phospho-Ser737) Antibody

Catalog No: #12427

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CFTR (Phospho-Ser737) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	CFTR (Phospho-Ser737) Antibody detects endogenous levels of CFTR only when phosphorylated at Ser737
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CFTR (Phospho-Ser737)
Target Name	CFTR
Modification	Phospho
Other Names	CFTR, ABC35, ABCC7, CF, DJ760C5.1, MRP7, CFTR/MRP, TNR-CFTR
Accession No.	Swiss-Prot#: P13569NCBI Gene ID: 1080
Uniprot	P13569
GeneID	1080;
Target Species	human
Calculated MW	168kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only