

Dihydrolipoyl dehydrogenase, mitochondrial

Catalog No: #AP82462



Package Size: #AP82462-1 50ug #AP82462-2 100ug #AP82462-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Dihydrolipoyl dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Dihydrolipoamide dehydrogenase
Accession No.	Uniprot ID:P49819
Uniprot	P49819
GeneID	403978;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DLD

Note: This product is for in vitro research use only